

LED Clear Skin Therapy



SOUND BITE: "Power up your clear skin treatments."

4-in-1 system for blemish-challenged skin combines four power-packed skin treatments into one tool to target blemishes, improve skin tone, minimize excess sebum, and infuse water-based treatments for clearer, healthier-looking skin. LED blue light, ultrasonic waves, and positive and negative ions gently and visibly improve the appearance of skin imperfections after every treatment.

WHAT DOES IT DO?

- **LED blue light** targets blemishes & clogged pores by eliminating acne-causing bacteria, improving clarity & minimizing the look of skin imperfections
- **Ultrasonic waves** stimulate blood circulation & improve the appearance of skin tone & texture
- **Positive ions** deep cleanse pores of excess oil & impurities
- **Negative ions** accelerate the absorption of serums & other water-based skincare products

DIRECTIONS

Apply water-based product, serum, or gel to targeted areas. Press power button on, then ION+ button, and move unit across face to cleanse pores. Press LED button, SONIC button, and ION- button, and move unit across targeted areas in circular motions, keeping skin moist. Do not use near the eyes or eyelids.

SKIN TYPES

all skin types

CONDITIONS

Blemishes, congested pores, enlarged pores, uneven tone, lines & wrinkles

CAUTION

Keep out of reach of children. Operate this unit only according to its intended use as described in the product manual. Do not use this unit if you have serious allergies, infected skin, or a heart condition. Do not immerse the unit in water. Keep the unit free of lint and hair.