Platinum Youth Restore & Contour

This lifting and contouring facial showcases the latest innovations from Platinum Pro Collection to visibly tighten, firm, counteract signs of environmental damage, and encourage recovery and renewal.

Platinum Pro Collection features concentrated formulas with Hya10+ Complex, a science-driven bio-cellular base containing 10 forms of hyaluronic acid boosted by powerful galactomyces ferment filtrate, designed to work with skin's own cellular systems for enhanced hydration, nourishment, and recovery.

Cleanse I	Perform first cleanse with PurErb Bare Beauty Cleansing Oil .
Analysis	Identify and confirm client skin concerns and assess skin barrier integrity.
Cleanse II	Dispense half pump Peptide Foam Cleanser into dish. Add water and whip with fan brush.
	Add 1 dropper AML 20 and mix.
	Apply to face, neck, and chest, avoiding orbital areas.
	Perform cleansing routine and remove.
Serum Infusion	Apply Iso-Rose Hydrating Mist generously.
	Apply 1 pump Beauty Reset Flash Firming Concentrate to face, neck, and chest.
Meridian	Using finger pads, press and hold Meridian points for 5 seconds each.
Massage	Apply 1 drop Derma Relief Serum to orbital areas. Apply as needed for slip and glide.
	Using finger pads or Mountain Shan Gua Sha Tool , press Meridian points, followed by a
	sweeping motion from center of the face and working outward.
Mask	Apply ProFerm Face Glow Mask to face.
	Apply ProFerm Neck Contour Mask to neck.
Finish	Choose Platinum Pro Targeted Serum:
	• To hydrate and brighten, apply 1 pump Youth Drops Daytime Firming Serum.
	• To calm and nourish, apply 1 pump Dream Drops Overnight Corrective Serum.
	Apply Bio Cell Rejuvenating Cream.
	Apply 1 pump Eye Wonder Youth Contour Eye Cream to orbital areas.
	Finish with Just Glow BB Cream SPF 50.

